5424

EN

Instructions for Use
Double lumen rectal balloon catheter, 10Fr

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Bruksanvisning
Dubbellumen rektal ballongkateter, 10 Fr

Intended use:

For use during urodynamics studies to measure the pressure in the rectum or vagina to determine the abdominal pressure. They are used in conjunction with bladder and filling catheters. They are used on a transient basis.

Environment of use: Hospital and sub-acute settings.

Patient population: Male and female patients with LUTS

Use with other devices:

Mediplus catheters use industry standard luers and connectors and can be used with any valve in line with the above intended use.

Contraindications:

1. Urinary Tract Infection
2. Indwelling Catheter
3. Known autonomic dysreflexia
4. Medication for Bladder dysfunction

Clinical Benefits:

To assist in the diagnosis of Lower Urinary Tract Symptoms (LUTS) such as bladder, urge, and stress urinary incontinence in conjunction with a urodynamic measurement system.

Performance Characteristics:

There are no specific performance claims made for Mediplus catheters. They are designed to perform as intended, e.g. to assist in the diagnosis of Lower Urinary Tract Symptoms, with minimal risks, comparable to similar devices. For the performance characteristic and product dimensions, refer to marketing brochures that are made available on the Mediplus website.

<https://www.mediplusuk.com/products>

Instructions:

- Attach red connector to pressure transducer on the urodynamics machine
- Open vent cap
- Hold catheter vertically with balloon end down and fill the catheter with fluid via the single red line to expel air
- When the fluid expels from the luer of the clear line, close cap
- Insert the catheter into the rectum or vagina
- Optional: evacuate fluid from balloon, using syringe in bottom of transducer, prior to insertion
- Refill balloon to desired volume after insertion

Warnings and Precautions

- This product is non-sterile and is intended for single use only

Caution: Federal Law (USA) restricts this device to sale by or on the order of a Physician

Device Disposal

The device should be disposed of as per the local healthcare policy and procedures. No special disposal treatment is required.

Serious incident notice:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer to the contact details included in this IFU and to the competent authority of the Member State where the serious incident occurred.

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Mediplus Ltd
Mediplus GmbH