STERILE EO

**5762, 5763, 5764, 5765,
5766, 5767, 5768, 5769**

EN

Instructions for Use

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Intended Use:

The 2-Way Silicone Foley Balloon Catheter is intended for use as urinary catheter to pass fluids to and from the urinary bladder. The Foley Catheters can be used for either Urethral or Suprapubic bladder drainage.

Environment of use: Hospital and sub-acute settings.

Patient Population: Patients requiring bladder irrigation and / or drainage.

Use with other devices: Mediplus urinary catheters use industry standard luers and connectors and can be used with any drainage bag or manual valve in line with the above intended use.

Contraindications:

Suprapubic use:

1. Carcinoma of the bladder
2. Anticoagulation and antiplatelet treatment
3. Abdominal wall sepsis
4. The presence of a subcutaneous vascular graft in the suprapubic region (e.g. a femoro-femoral crossover graft)

Urethral use:

1. Urethral injury or stricture
2. Recent urethral or bladder surgery, including but not limited to Urethroplasty, repair/reconstruction surgery, certain prostate procedures
3. Prior transurethral resection of the prostate with a large tissue defect
4. Significant symptoms of urinary obstruction prior to treatment
5. History of abdominoperineal resection for rectal cancer, rectal stenosis, or other major rectal pathology

Clinical Benefits:

Urinary catheters are tubes that are used to drain urine from the bladder for a range of medical reasons. An indwelling urinary catheter is kept in place by a balloon that is inflated once the catheter is in the bladder. A Suprapubic Catheter is used for drainage of the bladder through a needle puncture or a stab wound over the bladder.

Performance Characteristics:

There are no specific performance claims made for Mediplus Foley catheters and kits. They are designed to perform as intended, e.g. to treat urinary retention and incontinence, with minimal risks, comparable to similar devices. For the performance characteristic and product dimensions, refer to marketing brochures that are made available on the Mediplus website.

<https://www.mediplusuk.com/products>

Procedure:

Catheter insertion

1. Insert catheter according to relevant local and national policies. Ensure asepsis is maintained throughout insertion.
2. Insert catheter using aseptic non-touch technique (ANTT) technique. Hold catheter using clear inner packaging. Do not touch catheter directly.
3. Insert catheter into bladder. X-ray may be used to confirm the placement of the catheter within the bladder due to the radiopaque stripe and tip.
4. Once catheter is within the bladder, inflate catheter balloon, via Luer valve to the volume indicated on the catheter and packaging using pre-filled syringe of sterile water or 10% Aqueous Glycerine solution. **NB.** If there is any uncertainty regarding the position of the catheter do not inflate balloon as this may harm the patient.
5. Attach suitable sterile urine drainage device to catheter.
6. In order to minimize bladder neck and meatal pressure secure catheter or drainage bag using an appropriate fixation device, as per local policy.
7. Record procedure in patient record as per local policy. Include catheter traceability label in record.
8. Periodic observation of patient is required to ensure that the catheter is draining correctly and will allow for any potential problems to be identified and corrected.
9. Dispose of all used equipment according to local policy

Catheter Removal

1. Apply gloves to clean hands
2. Deflate catheter balloon using 10ml syringe. **NB.** If balloon volume in existing catheter is greater than 10ml, multiple aspirations will be required to ensure catheter balloon is completely empty
3. Gently withdraw catheter
4. Wipe meatus with gauze
5. Dispose of catheter and used equipment according to local policy.

Replacement Schedule

The schedule for replacing the catheter should be based on clinical symptoms for each individual patient or as often as recommended by a physician, but no longer than 29 days.

Warnings and precautions:

- Do not use ointments or lubricants made from a petroleum base
- In the event that the balloon will not deflate when the sterile water is removed, contact an adequately trained professional for assistance, as directed by the local healthcare policy and protocol
- In the event of a decrease in performance, patients are advised to contact a healthcare professional immediately
- Single Use Only
- Observe the expiry date on the sterilized peel pouch label prior to use
- Sterilised by ethylene oxide gas. Do **not** re-sterilize or reuse. Reuse of the catheter may lead to the risk of inflammation reaction and bacterial infection
- If using an accessory guide wire (not supplied) with the Open Tip catheters (5763, 5765, 5767 & 5769) then please observe all local healthcare policies and protocols and any relevant manufacturer's instructions
- Sterile unless package is opened or damaged
- Visually inspect that packaging is not damaged or unintentionally opened before use.

- Product should be discarded if the sterile packaging is damaged or unintentionally opened

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician

Device Disposal:

The device should be disposed of as per the local healthcare policy and procedures. No special disposal treatment is required.

Adverse events:

Suprapubic use

- Infection – including urinary tract infection and infection of the track site or wound.

Urethral use

- Urethral injury
- Skin breakdown
- Bladder stone
- Haematuria

Serious incident notice:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer to the contact details included in this IFU and to the competent authority of the Member State where the serious incident occurred

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Mediplus Ltd
Mediplus GmbH