STERILE EO

S-Cath™ System
5753, 5757, 5759,
5758, 5760, 5754, 5756

S-Cath™ System LONG
5759-L

EN

Instructions for Use
Suprapubic Foley Catheter Introducing Set

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CE **UK**
2797 **CA**₀₀₈₆

Product code 5758 is not for sale in the USA

Intended use:

The Suprapubic Catheter with Introducer Set is indicated for drainage of the bladder through a needle puncture or stab wound over the bladder. For use no greater than 29 days.

Environment of use: Hospital and sub-acute settings.

Patient population: Patients requiring bladder irrigation and / or drainage.

Use with other devices:

Mediplus urinary catheters use industry standard luers and connectors and can be used with any drainage bag or manual valve in line with the above intended use.

Contraindications to SPC Usage

- Carcinoma of the bladder
- Anticoagulation and antiplatelet treatment
- Abdominal wall sepsis
- The presence of a subcutaneous vascular graft in the suprapubic region (e.g. a femoro-femoral crossover graft)

Clinical Benefits:

Urinary catheters are tubes that are used to drain urine from the bladder for a range of medical reasons. An indwelling urinary catheter is kept in place by a balloon that is inflated once the catheter is in the bladder. A Suprapubic Catheter is used for drainage of the bladder through a needle puncture or a stab wound over the bladder.

Performance characteristics

There are no specific performance claims made for Mediplus Foley catheters and kits. They are designed to perform as intended, e.g. to treat urinary retention and incontinence, with minimal risks, comparable to similar devices. The Mediplus S-Cath system is safer than the competitors. For the performance characteristic and product dimensions, refer to marketing brochures that are made available on the Mediplus website.

<https://www.mediplusuk.com/products>

Preparation

Position the patient supine with a full bladder. Ensure bladder is palpable. The use of ultrasound is recommended to ensure that there is no entrapment of bowel between the bladder and abdomen.

Insertion Instructions

1. Connect the catheter to the urine drainage bag and prepare a syringe of 5-10ml sterile water to inflate the catheter balloon. If using 8Fr. Size (code 5758), prepare a syringe of 3ml sterile water.
2. Prepare the skin with antiseptic solution.
3. Administer local anaesthetic to the skin approximately 2 to 3cm above the Symphysis Pubis in the midline.
4. Make a small incision with the scalpel at this point.
5. Using the long needle and syringe inject a local anaesthetic along a track from the skin incision down to the bladder perpendicular to the skin. Continue to advance the needle until fluid from the bladder can be freely aspirated.
6. Remove the syringe of local anaesthetic from the needle and insert the guidewire through the needle, floppy end first, up to the first black mark.
7. Remove the needle being careful to hold the guidewire in place so it does not slip out of position. Feed the trocar over the guidewire ensuring that the guidewire is not pushed further into the bladder.
8. With a twisting action and steady gentle pressure insert the trocar into the bladder along the guidewire until the second black mark becomes visible. At this point the tip of the sheath will be in the bladder.
9. Remove the trocar and guidewire from the outer sheath. Some fluid will leak from the sheath. Occlude the outlet of the sheath with your thumb to avoid excessive release of fluid from the bladder.
10. Insert the Foley Catheter down the sheath to approximately the mid-point of the catheter length. If using 8Fr. Size (code 5758), remove internal guidewire from catheter before proceeding to next step.
11. Inflate the catheter balloon with 5-10ml* of sterile water. *If using 8Fr. Size (code 5758), inflate the catheter balloon with 3ml of sterile water.
12. Slide the outer sheath back along the catheter shaft until it is external to the abdomen. Pull the sheath apart and remove it from the catheter.

No external fixation is required.

The site of the catheter insertion can now be treated according to normal local procedures.

Please ensure disposal of needle, trocar and scalpel is in accordance with hospital Sharps policy.

Removal Instructions and Replacement Schedule

- Deflate the balloon completely and withdraw the catheter. Prepare the skin with antiseptic solution
- If catheterisation is to be discontinued, apply a dry dressing to the puncture site. If recatheterisation is necessary, follow normal local procedures
- The schedule for replacing the catheter should be based on clinical symptoms for each individual patient or as often as recommended by physician

Device Disposal

The device should be disposed of as per the local healthcare policy and procedures. No special disposal treatment is required.

Warnings and Precautions

- Do not use ointments or lubricants made from a petroleum base
- In the event that the balloon will not deflate when the sterile water is removed, contact an adequately trained professional for assistance, as directed by local healthcare policy and protocol
- In the event of a decrease in performance of the device, patients are advised to contact a healthcare professional immediately
- Single Use Only
- Observe the expiry date on the sterilized peel pouch label prior to use
- Sterilised with ethylene oxide gas. Do **not** re-sterilize or reuse. Reuse of the catheter may lead to the risk of inflammation reaction and bacterial infection
- Visually inspect that packaging is not damaged or unintentionally opened before use.
- Product should be discarded if the sterile packaging is damaged or unintentionally opened before use.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Adverse Events

- Haemorrhage – including haematuria and intra-abdominal bleeding
- Infection – including urinary tract infection and infection of the track site or wound
- Pain
- Injury to abdominal organs

Serious incident notice:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer to the contact details included in this IFU and to the competent authority of the Member State where the serious incident occurred.

Patented GB2410901 & GB2447163; International Patents applied for:
European Patent Application No. 05708317.2; US Patent Application No. 10/597,893.

Acknowledgement:
The Suprapubic Foley Catheter Introducing Set has been developed in conjunction with
Professor R.C.L. Feneley, of the Bristol Urological Institute, UK.

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Mediplus Ltd
Mediplus GmbH