



# SILICONE VESSEL ID LOOPS





# SILICONE VESSEL ID LOOPS

## INTENDED USE:

The “Fortune®” Silicone Vessel ID Loops are used in complicated surgical procedures by offering retraction, occlusion, and identification of arteries, nerves, veins, tendons or ureters.

## Fortune® Silicone Vessel ID Loops features:

- Medical grade silicone for superior biocompatibility.
- Soft and smooth material for better comfort
- Silicone radiopaque
- Non-adherent and non-absorbent to tissue
- Non-toxic and latex-free
- Sterile package
- Available in four sizes and four colors




## SPECIFICATIONS

REF. No.	Size.	Application
9005-0503R		Red Color
9005-0503B	Width: 0.5 mm Thickness: 0.3 mm Length: 450 mm	Blue Color
9005-0503Y		Yellow Color
9005-0503W		White Color

REF. No.	Size.	Application
9005-1309R		Red Color
9005-1309B	Width: 1.3 mm Thickness: 0.9 mm Length: 450 mm	Blue Color
9005-1309Y		Yellow Color
9005-1309W		White Color



9005-0503X



9005-1309X

GENERAL SURGICAL SERIES

# SILICONE VESSEL ID LOOPS

REF. No.	Size.	Application
9005-2510R		Red Color
9005-2510B	Width: 2.5 mm Thickness: 1.0 mm Length: 450 mm	Blue Color
9005-2510Y		Yellow Color
9005-2510W		White Color
9005-2510A		Red, Blue, Yellow, White Color

REF. No.	Size.	Application
9005-5013R		Red Color
9005-5013B	Width: 5.0 mm Thickness: 1.3 mm Length: 450 mm	Blue Color
9005-5013Y		Yellow Color
9005-5013W		White Color



9005-2510X



9005-5013X



**fortune  
medical**

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
Danshuei Dist, New Taipei City. 251 Taiwan  
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
www.fortunemed.com Mail:info@fortunemed.com



**Exclusive Distribution:**

Mediplus Ltd  
Unit 7, The Gateway Centre, Coronation Road  
Cressex Business Park, High Wycombe, Bucks, HP12 3SU  
Tel: 01494 551200  
Email: help@mediplusuk.com  
www.mediplusuk.com