



SILICONE HEMATURIA CATHETER

Your Right Choice for Quality Catheters



UROLOGICAL
SERIES



SILICONE HEMATURIA CATHETER

Fortune® Silicone 3-Way Hematuria Catheter Features:

- 100% medical grade silicone for superior biocompatibility.
- Transparent tube allows easy visual inspection and fluid observation.
- X-ray opaque line allows intubating location confirmation.
- A soft and uniformly inflated larger volume balloon makes the tube sit well against the bladder.
- 3-Way design with an attached plug on the irrigation funnel.
- Tapered Central funnel could be connected with irrigation bags.
- Five types are available for various applications:
 - A. type: Straight round tip
 - B. type: Straight whistle tip
 - C. type: Coude round tip
 - D. type: Dufour
 - E. type: Couvelaire
- Sterilized double packaging.


- A** Silicone Hematuria Catheter, Straight Round Tip
- B** Silicone Hematuria Catheter, Straight Whistle Tip
- C** Silicone Hematuria Catheter, Coude Round Tip
- D** Silicone Hematuria Catheter, Dufour
- E** Silicone Hematuria Catheter, Couvelaire



SPECIFICATIONS

Hematuria Catheter, Straight Round Tip:


REF. No.	Size	Balloon	Length	Description
1861-0318	18 FR	50-80cc		Straight round tip
1861-0320	20 FR		420 mm	X-ray opaque line
1861-0322	22 FR	80-100cc		Attached plug
1861-0324	24 FR			Catheter spigot



SILICONE HEMATURIA CATHETER


Hematuria Catheter, Straight Whistle Tip:

REF. No.	Size	Balloon	Length	Description
1862-0318	18 FR	50-80 cc	420 mm	Straight whistle tip X-ray opaque line Attached plug Catheter spigot
1862-0320	20 FR			
1862-0322	22 FR	80-100 cc		
1862-0324	24 FR			




Hematuria Catheter, Coude Round Tip:

REF. No.	Size	Balloon	Length	Description
1863-0318	18 FR	50-80 cc	420 mm	Coude round tip X-ray opaque line Attached plug Catheter spigot
1863-0320	20 FR			
1863-0322	22 FR	80-100 cc		
1863-0324	24 FR			




Hematuria Catheter, Dufour:

REF. No.	Size	Balloon	Length	Description
1864-0318	18 FR	50-80 cc	420 mm	Dufour X-ray opaque line Attached plug Catheter spigot
1864-0320	20 FR			
1864-0322	22 FR	80-100 cc		
1864-0324	24 FR			



Hematuria Catheter, Couvelaire:

REF. No.	Size	Balloon	Length	Description
1865-0318	18 FR	50-80 cc	420 mm	Couvelaire X-ray opaque line Attached plug Catheter spigot
1865-0320	20 FR			
1865-0322	22 FR	80-100 cc		
1865-0324	24 FR			



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com

Distributed by:

Mediplus Ltd., The Gateway Centre, Coronation Road, High Wycombe, HP12 3SU

T: +44 (0) 1494 551200
F: +44 (0) 1494 536333
help@mediplusuk.com



www.mediplusuk.com