




# Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics



# Mediplus Urodynamics Range


The Gold Standard in Fluid Filled Urodynamics

**For over 30 years Mediplus has manufactured a gold standard range of catheters and accessories designed to enhance safety and comfort during urodynamic procedures.**

**The Mediplus urodynamic range was developed in conjunction with clinicians to ensure optimum performance and ease of use.**

**Mediplus catheters offer a wide range of benefits including –**

- Colour coded for ease of identification
- Integrated extension sets to reduce the number of connections, thereby reducing inventory and lowering costs
- Ergonomically designed with winged luer locks for optimum grip and secure connection to transducers
- Specific lumen split reduces the size of the catheter whilst maintaining filling rates
- Smooth 'eyes' reduce discomfort when inserting and removing the catheter




# Mediplus Urodynamics Range


The Gold Standard in Fluid Filled Urodynamics

## Female Urodynamics

Many different configurations of a urodynamics "circuit" are used throughout the world. These illustrations show a typical UK/ European configuration however, whatever your requirements, Mediplus can supply a range of components to suit your needs.



## Male Urodynamics





# Mediplus Urodynamics Range

## The Gold Standard in Fluid Filled Urodynamics

### Urodynamics Catheters

Product Code	Product Description	Size	Length	Box Quantity
<b>Single Lumen Filling &amp; Cystometry Catheters</b>				
5100	Bladder Pressure Catheter	4.5 Fr./Ch.	200cm	25
5206	Bladder Filling Catheter	6 Fr./Ch.	40cm	10
5208	Bladder Filling Catheter	8 Fr./Ch.	40cm	10
5200	Bladder Filling Catheter	10 Fr./Ch.	40cm	10
5212	Bladder Filling Catheter	12 Fr./Ch.	40cm	10

Product Code	Product Description	Size	Filling	Line Length Pves	Double Lumen	Box Quantity
<b>Double Lumen Cystometry Catheters</b>						
5716	Cystometry Catheter	6 Fr./Ch.	15cm	130cm	65cm	10
5717	Cystometry Catheter (short)	6 Fr./Ch.	15cm	130cm	20cm	10
5736	Cystometry Catheter (standard)	6 Fr./Ch.	7cm	7cm	30cm	10
5702	Cystometry Catheter	8 Fr./Ch.	15cm	130cm	65cm	10
5703	Cystometry Catheter	8 Fr./Ch.	15cm	130cm	30cm	10
5732	Cystometry Catheter	8 Fr./Ch.	7cm	7cm	30cm	10
5712	Cystometry Catheter	10 Fr./Ch.	15cm	130cm	65cm	10





# Mediplus Urodynamics Range

## The Gold Standard in Fluid Filled Urodynamics

### Urodynamics Catheters

Product Code	Product Description	Size	Filling	Line Lengths Pves & Pura	Multi-lumen	Box Quantity
<b>Urethral Pressure Profile Catheters</b>						
5300	Single Lumen Catheter	8 Fr./Ch.	-	40cm	-	10
5724Y	Double Lumen Catheter with Y connector	6 Fr./Ch.	15cm	130cm	65cm	10
5704	Double Lumen Catheter	8 Fr./Ch.	15cm	15cm	65cm	10
5706	Double Lumen Catheter, extended lines	8 Fr./Ch.	130cm	130cm	65cm	10
5704Y	Double Lumen Catheter with Y connector	8 Fr./Ch.	15cm	130cm	65cm	10
5331	Triple Lumen Catheter	9 Fr./Ch.	15cm	130cm	50cm	10
5332	Triple Lumen Catheter (short), packed coiled	9 Fr./Ch.	15cm	15cm	50cm	10
5335	Triple Lumen Catheter (short), packed straight	9 Fr./Ch.	15cm	15cm	35cm	10




# Mediplus Urodynamics Range

## The Gold Standard in Fluid Filled Urodynamics

### Abdominal (Rectal) Pressure Catheters

Product Code	Product Description	Size	Balloon	Sterile	Length	Box Quantity
<b>Abdominal (Rectal Pressure Catheters)</b>						
5145	Rectal Catheter, Sterile	4.5 Fr./Ch.	-	✓	200cm	100
5400	Rectal Catheter, Sterile	4.5 Fr./Ch.	Flap valve	✓	200cm	25
5405	Rectal Catheter	4.5 Fr./Ch.	Flap valve	✗	200cm	25
5410	Rectal Balloon Catheter, Sterile	4.5 Fr./Ch.	With slit	✓	200cm	25
5415	Rectal Catheter	4.5 Fr./Ch.	With slit	✗	200cm	25
5435	Rectal Balloon Catheter, non-sterile	4.5 Fr./Ch.	(small) with slit	✗	200cm	25
5415i	Rectal Balloon Catheter with introducer, non-sterile	4.5 Fr./Ch.	With slit	✗	200cm	25
5457	Rectal Balloon Catheter, non-sterile	10 Fr./Ch.	With slit	✗	38.5cm	10
5424	Double Lumen Rectal Balloon Catheter, non-sterile*	10 Fr./Ch.	✓	✓	38.5cm	10

\*Cystometry lines provide additional 130cm length



# Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics

## Extension Sets

Product Code	Product Description	Tubing i.d./o.d.	Length	Box Quantity
9003	Male/female luer lock	1.5mm / 3.10mm	5cm	50
9015	Male/female luer lock	1.5mm / 3.10mm	15cm	50
9025	Male/female luer lock	1.5mm / 3.10mm	25cm	50
9050	Male/female luer lock	1.5mm / 3.10mm	50cm	50
9100	Male/female luer lock	1.5mm / 3.10mm	100cm	50
9111	Double extension set (Siamese tube)	1.5mm / 3.10mm	100cm	20
9115	Double extension set (Siamese tube)	1.5mm / 3.10mm	150cm	20
9200	Male/female luer lock	1.5mm / 3.10mm	200cm	50
5183	Extension Set with one-way valve	3.0mm / 4.10mm	20cm	20

## Accessories

Product Code	Product Description	Tubing i.d./o.d.	Length	Box Quantity
5182	SetGuard	3.5mm / 5.50mm	30cm	25



Double extension set (Siamese tube)



Extension Set with one-way valve



SetGuard



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.  
+44 (0) 1494 551200  
help@mediplusuk.com  
The Gateway Centre, Coronation Road,  
High Wycombe HP12 3SU  
www.mediplusuk.com

**INVESTORS IN PEOPLE™**  
We invest in people Silver

