



Mediplus POPY™ Pessary

A conservative, cost-effective solution for pelvic organ prolapse or procidentia



Mediplus POPY™ Pessary

A conservative, cost-effective solution for pelvic organ prolapse or procidentia

The Mediplus POPY™ was the world's first silicone shelf pessary. This innovative pessary offers a wide range of benefits including:

Conservative management – POPY™ is a conservative option for managing pelvic organ prolapse

Ease of use – high-quality medical grade silicone ensures that the POPY™ is soft and malleable for ease of fitting and removal

Reduced risk – the softness of the POPY™ pessary reduces the risk of erosion

Low-cost – compared to surgical intervention

Innovative design – the POPY™ range has been developed with patients' needs at the forefront of its design




“POPY™ is a far better choice for the patient as it is not rigid at all but rather flexible, therefore insertion and removal is not painful for my patients”

Claudine Domoney, Lead Urogynae Consultant,
Chelsea & Westminster Hospital


Mediplus POPY™ Pessary

A conservative, cost-effective solution for pelvic organ prolapse or procidentia


Insertion technique




1. Prolapse



2. Initial insertion



3. Pessary inserted further



4. Pessary in situ



Ordering information

Product Code	Product Description	Size	Width (mm)	Box Quantity
POPY57	POPY™ Silicone Shelf Pessary	6	57	1
POPY60	POPY™ Silicone Shelf Pessary	7	60	1
POPY64	POPY™ Silicone Shelf Pessary	8	64	1
POPY67	POPY™ Silicone Shelf Pessary	9	67	1
POPY70	POPY™ Silicone Shelf Pessary	10	70	1
POPY73	POPY™ Silicone Shelf Pessary	11	73	1
POPY76	POPY™ Silicone Shelf Pessary	12	76	1
POPY79	POPY™ Silicone Shelf Pessary	13	79	1
POPY83	POPY™ Silicone Shelf Pessary	14	83	1
POPY86	POPY™ Silicone Shelf Pessary	15	86	1
POPY89	POPY™ Silicone Shelf Pessary	16	89	1
POPY57-S	POPY™ Silicone Shelf Pessary - short stem	6	57	1
POPY60-S	POPY™ Silicone Shelf Pessary - short stem	7	60	1
POPY64-S	POPY™ Silicone Shelf Pessary - short stem	8	64	1
POPY67-S	POPY™ Silicone Shelf Pessary - short stem	9	67	1
POPY70-S	POPY™ Silicone Shelf Pessary - short stem	10	70	1
POPY73-S	POPY™ Silicone Shelf Pessary - short stem	11	73	1
POPY76-S	POPY™ Silicone Shelf Pessary - short stem	12	76	1
POPY79-S	POPY™ Silicone Shelf Pessary - short stem	13	79	1
POPY83-S	POPY™ Silicone Shelf Pessary - short stem	14	83	1
POPY86-S	POPY™ Silicone Shelf Pessary - short stem	15	86	1
POPY89-S	POPY™ Silicone Shelf Pessary - short stem	16	89	1

Clinical data

- ¹ One year after treatment women with POP report similar improvement when treated with a pessary versus surgical correction. Prospective evaluation of outcome of vaginal pessaries versus surgery in women with symptomatic pelvic organ prolapse. Z. Abdool, R.Thakar, A.H Sultan, R.S Oliver. Int Urogynecol J (2011) 22: 273-278
- ² POP management with pessaries is shown to have a significant improvement in enjoyment of life and interference with physical activity. Management of pelvic organ prolapse by vaginal pessaries - one year follow up. Oliver R, Lim J H, Thakar R, Sultan A H. Mayday University Hospital, Croydon, Surrey, UK.
- ³ Pessaries have been shown to be as effective as surgery and therefore should be offered as an alternative option to surgery. Are vaginal pessaries as effective as surgery in symptomatic pelvic organ prolapse? Z.Abdool, R. Thakar, A.Sultan. Int Urogynecol J (2006) 17 (Suppl.2)
- ⁴ Vaginal pessaries are a cost effective way of alleviating symptoms of pelvic organ prolapse. Guidelines for the use of support pessaries in the management of pelvic organ prolapse University of South Australia (2010)
- ⁵ Vaginal pessary use is a non-invasive treatment option for patients with POP who wish to postpone surgery or avoid surgical risks altogether. Pelvic Organ Prolapse - Pessary Treatment J.A Schulz, E. Kwon. American Journal of Obstetrics and Gynecology 2005;193 (1): 103-13
- ⁶ Pessaries provide mechanical support to the prolapsed organs, thus relieving symptoms and improving overall body image for the patient. Pessary use and impact on quality of life and body image. Patel M.S, Mellen C, O'Sullivan D. M, Lasala C.A. Female Pelvic Medical Reconstructive Surgery (2011) Nov;17 (6): 298-301



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.

+44 (0) 1494 551200

help@mediplusuk.com

The Gateway Centre, Coronation Road,
High Wycombe HP12 3SU

www.mediplusuk.com

INVESTORS IN PEOPLE™
We invest in people Silver

