


USING TIVA? **PERUSE** BEFORE YOU *INFUSE*

P: PROPOFOL

Propofol in the syringe & “Propofol” on the pump?

Is this the intended model?




E: EMPTY AMPOULES

Empty remifentanyl ampoule?

1mg or 2mg?

Correct dilution?


Correct propofol concentration?



R: REMIFENTANIL

Remifentanyl in the syringe & “Remifentanyl” on the pump?


Is this the intended model?



U: USER SETTINGS

Are the patient characteristics the same?

Are they THIS patient?



S: SETUP / GIVING SET

New, appropriate giving set?

Primed, unclamped & connected at BOTH ends?


Cannula visible & drip running?



E: ELECTRICITY

Do the pumps have adequate charge to begin?

Are there power cables ready for use?



TWO-PERSON POSITIVE CONFIRMATORY CHECK, BEFORE STARTING THE INFUSIONS

Version 1.4. R. Fleming 2022
QR code links to Association of Anaesthetists “Education Shot”
Pump images with permission © Becton, Dickinson & company

