

A collection of Mediplus peripheral IV connectors. One connector is shown in the upper left, attached to a clear tube with a green handle and a red port. Below it, several other connectors in various colors (clear, blue, green, red) are shown, some with red ports. One clear connector has 'GRAVITY IV' printed on it. The background is white with a green horizontal band across the middle.

Mediplus Peripheral IV Connectors

Safer, more efficient drug delivery

Mediplus Peripheral IV Connectors

Safer, more efficient drug delivery

Mediplus peripheral IV connectors are designed to ensure optimal drug delivery for enhanced safety in ICU, HDU, orthopaedic and obstetric departments. These high-performance connectors offer a wide range of benefits including:

Enables rapid emergency response – the high-flow anti-reflux valve (Coventry Connector and High Flow Y Connector) offers a uniquely fast flow rate for rapid IV fluid delivery and prevents back flow into the IV line

Greater control – the integrated anti-siphon valve prevents free flow and back flow, reducing the risk of over or under dosing during infusion


Comfort & convenience – the removable bridge (Coventry and Chalfont Connectors) helps prevent the device becoming dislodged and keeps it compact for greater patient comfort and manoeuvrability

Ergonomic design – the J-loop design (Coventry and Chalfont Connectors) prevents kinking and occlusion of tubing for unimpeded flow

Wide range – single, double or triple infusion options available

“I have been using the Coventry Connectors for more than 10 years now and I find them really useful. The multiple ports make it easy to run multiple infusions and fluids together. They are kink resistant and the rate of fluid flow is really good”

Payal Kajekar,
Luton & Dunstable Hospital



Mediplus Peripheral IV Connectors

Safer, more efficient drug delivery

Key features



J-loop to prevent kinking/occlusion of the tubing



Removable bridge offers more flexibility when moving patients



Anti-siphon valves to prevent back flow and free flow



Wide bore gravity line to allow high flow of IV fluid



Anti-reflux valve to prevent back flow



Reduced tube length, improving high flow rate


Peripheral IV connector range:



Chalfont Connector



Coventry Connector



High Flow Y Connector

Ordering information

Product Code	Product Name	Product Description	Box Quantity
6507	Chalfont	Double lumen extension set: 1 x anti-siphon, 1 x high-flow back check valve (13cm/5")	50
6510	Coventry	Triple lumen extension set: 2 x anti-siphon, 1 x high-flow back check valve (17cm/7")	50
6517	High Flow Y	Double lumen extension set: 2 x high-flow back check valve (13cm/5")	50

Clinical data

¹ Safe Anaesthesia Liaison Group (SALG), Guaranteeing Drug Delivery in Total Intravenous Anaesthesia, October 2009



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.
+44 (0) 1494 551200
help@mediplusuk.com
The Gateway Centre, Coronation Road,
High Wycombe HP12 3SU
www.mediplusuk.com

INVESTORS IN PEOPLE®
We invest in people Gold

