

A clear plastic capnomask with a blue strap and a green tube connected to a blue tube. The device is shown in a 3/4 view, highlighting its ergonomic shape and the connection points. The background is white with a green horizontal band across the middle.

# Mediplus Capnomask™

Fast, reliable EtCO<sub>2</sub> monitoring



# Mediplus Capnomask™

Fast, reliable EtCO<sub>2</sub> monitoring

## The high-performance Capnomask™ offers a wide range of benefits including:

**Enhanced patient safety** – provides accurate monitoring to help prevent hypoxia post-surgery

**Empowers fast clinical response** – the unique, front dual connector features fully integrated oxygen and carbon dioxide lines, enabling fast, accurate detection of end-tidal CO<sub>2</sub> and supporting swift clinical response

**Supports best practice** – the Association of Anaesthetists (AoA) recommends that continuous capnography should be used in all anaesthetised patients<sup>1</sup>

**Reduced risk of infection** – integrated 0.22 micron filter reduces the risk of equipment cross-contamination and patient infection

**Greater patient comfort** – ergonomically moulded in soft, high-quality PVC for a stable, comfortable fit

**Helps support faster recovery** – impregnated with a pleasant, subtle vanilla scent to help reduce post-operative nausea and vomiting

**Wide range of sizes** – adult, paediatric and ophthalmic versions available



**“The Capnomask™ has a superior attachment of the O<sub>2</sub> tubing to the mask and the mask quality is also superior and offers a great fit. I prefer it to any other mask on the market”**

Dr Emad Aziz - Guy's and St Thomas'  
NHS Foundation Trust


# Mediplus Capnomask™

Fast, reliable EtCO<sub>2</sub> monitoring

## Key features



Unique dual adapter enables accurate detection of EtCO<sub>2</sub>



Consistent monitoring at 10l/min



Unique vanilla scent supports recovery from PONV



Standard connector for oxygen outlet



Microstream™ connector option available for preferred CO<sub>2</sub> monitor compatibility also including integrated filter



Standard Luer connector with integrated filter protecting both patient and monitor from cross-contamination

## Ordering information

Product Code	Product Description	Filter	Box Quantity
A202Mx	Capnomask™ Adult Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202CMx	Capnomask™ Ophthalmic Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202PMx	Capnomask™ Paediatric Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202Fx	Capnomask™ Adult Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202CFx	Capnomask™ Ophthalmic Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202PFx	Capnomask™ Paediatric Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202Mx-D	Capnomask™ DISS Connector Adult Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202CMx-D	Capnomask™ DISS Connector Ophthalmic Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202PMx-D	Capnomask™ DISS Connector Paediatric Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/male connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202Fx-D	Capnomask™ DISS Connector Adult Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202CFx-D	Capnomask™ DISS Connector Ophthalmic Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202PFx-D	Capnomask™ DISS Connector Paediatric Oxygen Mask + CO <sub>2</sub> sampling line 300cm with male/female connector + O <sub>2</sub> tube 200cm	0.22μ	50
A202Ax-M	Capnomask™ Adult Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Cx-M	Capnomask™ Ophthalmic Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Px-M	Capnomask™ Paediatric Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Ax-D-M	Capnomask™ DISS Adult Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Cx-D-M	Capnomask™ DISS Ophthalmic Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Px-D-M	Capnomask™ DISS Paediatric Oxygen Mask + CO <sub>2</sub> Microstream™ sampling line 300cm + O <sub>2</sub> tube 200cm	0.22μ	50
A202Ac	Capnomask™ Adult Oxygen Mask	-	50

## Clinical data

<sup>1</sup> "Capnography monitoring is essential at all times in patients with tracheal tubes, supraglottic airway devices and those who are sedated to a level unresponsive to verbal commands." Updated statement from the Association of Anaesthetists (AoA) Guidelines. May 2021. <https://anaesthetists.org/Portals/0/PDFs/Guidelines%20PDFs/Recommendations%20for%20standards%20of%20monitoring%20during%20anaesthesia%20and%20recovery%202021.pdf?ver=2021-05-26-141701-007>  
Updated May 2021.



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.

+44 (0) 1494 551200

help@mediplusuk.com

The Gateway Centre, Coronation Road,  
High Wycombe HP12 3SU

www.mediplusuk.com

**INVESTORS IN PEOPLE®**  
We invest in people Gold

