

A clear plastic capnomask with a blue strap and a green tube connected to a blue tube. The device is shown in a side profile, highlighting its ergonomic shape and the connection points for the tubes.

Mediplus Capnomask™

Fast, reliable EtCO₂ monitoring

40
Years of innovative patient care

Mediplus Capnomask™

Fast, reliable EtCO₂ monitoring

The high-performance Capnomask™ offers a wide range of benefits including:

Enhanced patient safety – provides accurate monitoring to help prevent hypoxia post-surgery

Empowers fast clinical response – the unique, front dual connector features fully integrated oxygen and carbon dioxide lines, enabling fast, accurate detection of end-tidal CO₂ and supporting swift clinical response

Supports best practice – the Association of Anaesthetists (AoA) recommends that continuous capnography should be used in all anaesthetised patients¹

Reduced risk of infection – integrated 0.22 micron filter reduces the risk of equipment cross-contamination and patient infection

Greater patient comfort – ergonomically moulded in soft, high-quality PVC for a stable, comfortable fit

Helps support faster recovery – impregnated with a pleasant, subtle vanilla scent to help reduce post-operative nausea and vomiting

Wide range of sizes – adult, paediatric and ophthalmic versions available



“The Capnomask™ has a superior attachment of the O₂ tubing to the mask and the mask quality is also superior and offers a great fit. I prefer it to any other mask on the market”

Dr Emad Aziz - Guy's and St Thomas'
NHS Foundation Trust


Mediplus Capnomask™

Fast, reliable EtCO₂ monitoring

Key features



Unique dual adapter enables accurate detection of EtCO₂



Consistent monitoring at 10l/min



Unique vanilla scent supports recovery from PONV



Standard connector for oxygen outlet



Microstream™ connector option available for preferred CO₂ monitor compatibility also including integrated filter



Standard Luer connector with integrated filter protecting both patient and monitor from cross-contamination

Ordering information

Product Code	Product Description	Filter	Box Quantity
A202Mx	Capnomask™ Adult oxygen mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202CMx	Capnomask™ Ophthalmic oxygen mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202PMx	Capnomask™ Paediatric Oxygen mask + CO ₂ sampling line 300cm with male/male connector + O ₂ tube 200cm	0.22µ	50
A202Ax-M	Capnomask™ Adult Oxygen Mask + CO ₂ Microstream™ sampling line 300cm + O ₂ tube 200cm	0.22µ	50
A202Cx-M	Capnomask™ Ophthalmic Oxygen Mask + CO ₂ Microstream™ sampling line 300cm + O ₂ tube 200cm	0.22µ	50
A202Px-M	Capnomask™ Paediatric Oxygen Mask + CO ₂ Microstream™ sampling line 300cm + O ₂ tube 200cm	0.22µ	50
A202Ac	Capnomask™, mask only, adult	-	50

Clinical data

¹ "Capnography monitoring is essential at all times in patients with tracheal tubes, supraglottic airway devices and those who are sedated to a level unresponsive to verbal commands." Updated statement from the Association of Anaesthetists (AoA) Guidelines. May 2021. <https://anaesthetists.org/Portals/0/PDFs/Guidelines%20PDFs/Recommendations%20for%20standards%20of%20monitoring%20during%20anaesthesia%20and%20recovery%202021.pdf?ver=2021-05-26-141701-007>
Updated May 2021.



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.
+44 (0) 1494 551200
help@mediplusuk.com
The Gateway Centre, Coronation Road,
High Wycombe HP12 3SU
www.mediplusuk.com

INVESTORS IN PEOPLE
We invest in people Gold

