

Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics

40

Years of innovative patient care

Mediplus Urodynamics Range


The Gold Standard in Fluid Filled Urodynamics

For over 30 years Mediplus has manufactured a gold standard range of catheters and accessories designed to enhance safety and comfort during urodynamic procedures.

The Mediplus urodynamic range was developed in conjunction with clinicians to ensure optimum performance and ease of use.

Mediplus catheters offer a wide range of benefits including –

- Colour coded for ease of identification
- Integrated extension sets to reduce the number of connections, thereby reducing inventory and lowering costs
- Ergonomically designed with winged luer locks for optimum grip and secure connection to transducers
- Specific lumen split reduces the size of the catheter whilst maintaining filling rates
- Smooth 'eyes' reduce discomfort when inserting and removing the catheter




Mediplus Urodynamics Range


The Gold Standard in Fluid Filled Urodynamics

Female Urodynamics

Many different configurations of a urodynamics “circuit” are used throughout the world. These illustrations show a typical UK/ European configuration however, whatever your requirements, Mediplus can supply a range of components to suit your needs.



Male Urodynamics




Mediplus Urodynamics Range


The Gold Standard in Fluid Filled Urodynamics

Urodynamics Catheters

Product Code	Product Description	Size	Length	Box Quantity
Single Lumen Filling & Cystometry Catheters				
5100	Bladder Pressure Catheter	4.5 Fr./Ch.	200cm	25
5206	Bladder Filling Catheter	6 Fr./Ch.	40cm	10
5208	Bladder Filling Catheter	8 Fr./Ch.	40cm	10
5200	Bladder Filling Catheter	10 Fr./Ch.	40cm	10
5212	Bladder Filling Catheter	12 Fr./Ch.	40cm	10



Product Code	Product Description	Size	Filling	Line Length Pves	Double Lumen	Box Quantity
Double Lumen Cystometry Catheters						
5716	Cystometry Catheter	6 Fr./Ch.	15cm	130cm	65cm	10
5717	Cystometry Catheter (short)	6 Fr./Ch.	15cm	130cm	20cm	10
5736	Cystometry Catheter (standard)	6 Fr./Ch.	7cm	7cm	30cm	10
5702	Cystometry Catheter	8 Fr./Ch.	15cm	130cm	65cm	10
5703	Cystometry Catheter	8 Fr./Ch.	15cm	130cm	30cm	10
5732	Cystometry Catheter	8 Fr./Ch.	7cm	7cm	30cm	10
5712	Cystometry Catheter	10 Fr./Ch.	15cm	130cm	65cm	10





Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics

Urodynamics Catheters

Product Code	Product Description	Size	Filling	Line Lengths Pves & Pura	Multi-lumen	Box Quantity
Urethral Pressure Profile Catheters						
5300	Single Lumen Catheter	8 Fr./Ch.	-	40cm	-	10
5724Y	Double Lumen Catheter with Y connector	6 Fr./Ch.	15cm	130cm	65cm	10
5704	Double Lumen Catheter	8 Fr./Ch.	15cm	15cm	65cm	10
5706	Double Lumen Catheter, extended lines	8 Fr./Ch.	130cm	130cm	65cm	10
5704Y	Double Lumen Catheter with Y connector	8 Fr./Ch.	15cm	130cm	65cm	10
5331	Triple Lumen Catheter	9 Fr./Ch.	15cm	130cm	50cm	10
5332	Triple Lumen Catheter (short), packed coiled	9 Fr./Ch.	15cm	15cm	50cm	10
5335	Triple Lumen Catheter (short), packed straight	9 Fr./Ch.	15cm	15cm	35cm	10




Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics

Abdominal (Rectal) Pressure Catheters

Product Code	Product Description	Size	Balloon	Sterile	Length	Box Quantity
Abdominal (Rectal Pressure Catheters)						
5145	Rectal Catheter, Sterile	4.5 Fr./Ch.	-	✓	200cm	100
5400	Rectal Catheter, Sterile	4.5 Fr./Ch.	Flap valve	✓	200cm	25
5405	Rectal Catheter	4.5 Fr./Ch.	Flap valve	✗	200cm	25
5410	Rectal Balloon Catheter, Sterile	4.5 Fr./Ch.	With slit	✓	200cm	25
5415	Rectal Catheter	4.5 Fr./Ch.	With slit	✗	200cm	25
5435	Rectal Balloon Catheter, non-sterile	4.5 Fr./Ch.	(small) with slit	✗	200cm	25
5415i	Rectal Balloon Catheter with introducer, non-sterile	4.5 Fr./Ch.	With slit	✗	200cm	25
5457	Rectal Balloon Catheter, non-sterile	10 Fr./Ch.	With slit	✗	38.5cm	10
5424	Double Lumen Rectal Balloon Catheter, non-sterile*	10 Fr./Ch.	✓	✗	38.5cm	10

*Cystometry lines provide additional 130cm length



Mediplus Urodynamics Range

The Gold Standard in Fluid Filled Urodynamics

Extension Sets

Product Code	Product Description	Tubing i.d./o.d.	Length	Box Quantity
9003	Male/female luer lock	1.5mm / 3.10mm	5cm	50
9015	Male/female luer lock	1.5mm / 3.10mm	15cm	50
9025	Male/female luer lock	1.5mm / 3.10mm	25cm	50
9050	Male/female luer lock	1.5mm / 3.10mm	50cm	50
9100	Male/female luer lock	1.5mm / 3.10mm	100cm	50
9111	Double extension set (Siamese tube)	1.5mm / 3.10mm	100cm	20
9115	Double extension set (Siamese tube)	1.5mm / 3.10mm	150cm	20
9200	Male/female luer lock	1.5mm / 3.10mm	200cm	50
5183	Extension Set with one-way valve	3.0mm / 4.10mm	20cm	20

Accessories

Product Code	Product Description	Tubing i.d./o.d.	Length	Box Quantity
5182	SetGuard	3.5mm / 5.50mm	30cm	25



Double extension set (Siamese tube)



Extension Set with one-way valve



SetGuard



Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



anaesthesia



gynaecology



operating room



urology

Mediplus Ltd.
+44 (0) 1494 551200
help@mediplusuk.com
The Gateway Centre, Coronation Road,
High Wycombe HP12 3SU
www.mediplusuk.com

INVESTORS IN PEOPLE®
We invest in people Gold

